


Parts Manual

Oxford® Ascend



| NO | PART NUMBER | DESCRIPTION | QTY |
|----|-------------|--|-----|
| 1 | 0Y0066 | 100mm rear castor (braked) | 2 |
| 2 | 0Y0263 | Straight-line steering device | 1 |
| 3 | 0Y0073 | 100mm front castor inc. load bearing washer and tightening wrench (wrench not shown) | 2 |
| 4 | 0Y0206 | Mast extrusion including labels * | 1 |
| | 0Y0308 | Mast extrusion including labels ** | 1 |
| 5 | 0Y0207 | Boom extrusion - Cow horn style (including boom end cap) * | 1 |
| | 0Y0388 | Boom extrusion - Cow horn style (including boom end cap) ** | 1 |
| 6 | 0Y0208 | Push handle assembly (inc. bolts & fixings) | 1 |
| 7 | 0Y0103 | Casted boom to mast pivot joint (inc. connection bolts/fixings/spacers) | 1 |
| 8 | 0Y0092 | Base assembly (inc. labels, base casting and bottom plate, not internal parts) | 1 |
| 9 | 0Y0209 | Left leg extrusion (inc. fixings) | 1 |
| 10 | 0Y0210 | Right leg extrusion (inc. fixings) | 1 |
| 11 | 39001370 | Control box (Smart™ Monitor) | 1 |
| | 0Y0211 | Control box (Standard) | 1 |
| 12 | 11012180 | Battery pack | 1 |
| 13 | 39001309 | User handset (Smart™ Monitor only) | 1 |
| | 0Y0075 | User handset (Standard control box) | 1 |
| - | 39001179 | Calibration handset (Smart™ Monitor only) (not shown) | 1 |
| - | 39001178 | Actuator reset handset (Smart™ Monitor only) (not shown) | 1 |
| 14 | 39001316 | Smart™ Monitor retrofit kit (inc. control box, battery and user handset) | 1 |
| 15 | 39001183 | Charger and stand assembly (with UK charger lead) | 1 |
| 16 | 11012183 | Charger lead (UK) | 1 |
| 17 | 11012184 | Charger lead (EU) | 1 |
| 18 | 0Y0212 | Actuator (inc. location pins, sleeves, washers, nuts, labels) * [OBSOLETE] | 1 |
| | 0Y0306 | Actuator (inc. location pins, sleeves, washers, nuts, labels) ** | 1 |
| 19 | 0Y0069 | Battery and controller bracket (inc. fixings) | 1 |
| 20 | 0Y0093 | Leg opening actuator - 35mm stroke (inc. bolts and fixings) | 1 |
| 21 | 0Y0126 | Pivot caps | 2 |
| 22 | 0Y0102 | Fixings kit (all bolts, screws and fixings) | 1 |
| 23 | 0Y0077 | Mast to base electrical connection kit (inc. wires) | 1 |
| 24 | 0Y0078 | Foot push pad assembly complete | 2 |
| 25 | 0Y0112 | Mast locking knob assembly (inc. plastic base insert) | 1 |
| 26 | 0Y0632 | Knee pad assembly (inc. all bolts, fixings and leg strap) | 1 |
| 27 | 0Y0101 | Foot plate assembly | 1 |
| 28 | 0Y0129 | Base to actuator electrical connection kit (inc. wires) | 1 |
| 29 | 0Y0108 | Boom to mast pivot plastic caps | 2 |
| 30 | 0Y0076 | Mast side caps | 2 |
| 31 | SA/2696 | Sling retainer clips (Actuator/Boom) | 2 |
| - | 0Y1061 | Sling retainer clips (Cow Horn) - D ring (not shown) | 10 |
| - | 0Y0687 | Sling retainers - Oblong (not shown) | 10 |
| - | 0Y0213 | Decal set (no serial number labels) (not shown) | 1 |
| - | 0Y0070 | Handset clip (not shown) | 2 |
| - | 0Y0128 | Castor tightening wrench (not shown) | 1 |

* Last serial no. supplied: **0703L1146**. Please check with Customer Services prior to ordering.

** Please check with Customer Services prior to ordering.

